

HYGISEAT self-cleaning toilet

Ref. SUP1500



The HYGISEAT toilet washes, disinfects, dries the seat after each use and the bowl flushes automatically.

Works with a simple approach of the hand above the appliance, without physical contact.

Rotation of the seat cycle : water pressure wash mixed with Germibloc disinfectant, drying the seat after each use.

Cycle time: 12 seconds

Total safety

- Operates on 12 v DC with external 230/12V transformer.
- The seats stops rotating immediately in case of manual interference.
- Manual flush command button during power cuts.
- EC. 3 years pieces warranty.

Installation on GEBERIT/HYGISEAT mounting frame with concealed tank

Concealed fixing



OPTIONS :

- **SUP1065** : double detection per approach of the hand or cycle auto-start after user departure.
- **SUP1080** : spray of blue cleaning foam product into the WC on users arrival (SUP1065 included)
- **NEW / SUP1070** : efficient odors suction system. Rejection of contaminated air directly into the drain (SUP1065 included)
- **SUP1095** : WC in colours, using RAL colour chart

GERMIBLOC™ disinfectant product

Biodegradable dry product presented in a 72 g plastic bottle, easy to manipulate and store.

Germibloc™ disinfectant product replacement when the toilet beeps 3 times at the end of cycle, every 1000 uses. Norm NFEN 1275 Fungicide - Norm NFEN 1040 Bactericide

